



Bellerophon Biotec Pvt. Ltd.

Bellerophon Biotec Pvt. Ltd.



Sensi Drive - **VALUE LIFE WHILE YOU DRIVE**

A unique, self-testing, disposable, simple, inexpensive and laboratory calibrated alcohol breathalyzer, which detects the presence of alcohol in the body and gives an accurate indication of whether the blood alcohol content has exceeded a specific limit.

SAFE DRIVE PRODUCT CATALOGUE

Users of SensiDrive-

DRIVERS: All drivers need to check for their blood alcohol level if they are drunk. This shall save them from penal action as well as protect their life and vehicle.

PARENTS: Checking for teenage alcohol consumption.

DRIVING SCHOOLS: Checking students and/or instructors for the presence of alcohol.

MOTORSPORT: Promoters/organizers could check for the blood alcohol content amongst competitors.

BUS/TRUCK COMPANIES: Checking drivers for the presence of alcohol before they depart.

COMPANIES: Checking for alcohol consumption by workers, especially if workers are using vehicles, machinery, hazardous goods, power tools, working in hazardous areas, working at height, working with electricity/explosives/flammable goods or any other safety risk.

SCHOOLS/COLLEGES/UNIVERSITIES: Checking students for the presence of alcohol, possibly at school formals, dances, school camps and sporting events.

Advantages of SensiDrive-

- ✓ Quick and Accurate
- ✓ Low cost
- ✓ Easy to use
- ✓ Hygienic & reliable
- ✓ Self-assessment by driver
- ✓ Portable and Disposable
- ✓ Requires no batteries or power source
- ✓ Approved and certified for accuracy internationally (FDA, USA; STANDARDS, AUSTRALIA, FRENCH NF & GERMAN TUV).
- ✓ Saves from penal procedures and accidents
- ✓ Shelf life: 2 years

Digital Breath Analyzer V/s Disposable Breath Analyzer

A COMPARISON CHART


S.No.	Digital Breath Analyzer	Disposable Breath Analyzer
1.	High risk of infections such as T.B, Influenza and other respiratory tract infections	Zero infection Risk
2.	Quite heavy	Easy to carry
3.	Calibration required	Self calibrated
4.	Not environment friendly	Eco friendly
5.	Mouthpieces required	No Mouthpieces required
6.	Costly	Affordable alcohol detection device
7.	Sensor is required that needs regular maintenance	No sensor required
8.	Battery needs to be changed	No battery needed
9.	Cannot be used by large number of people	Can be used by masses
10.	Chances of false positive test exist	No chance of false positive test

How to Use SensiDrive-?

Test yourself within 2 minutes


STEP 1 : Snap-off End Caps

Break off the snap-off end caps at both ends of the tester and shake out the harmless white crystals inside the tube.



STEP 2 : Inflate the Bag

Inflate the bag by blowing in one steady breath into the neck of the bag. The valve in the neck of the bag automatically closes after inflation.



STEP 3 : Insert the White Cap

Insert the white cap end of the tester into the neck of the bag.




STEP 4 : Squeeze the Bag

Squeeze the bag evenly, expelling the air through the tester.



STEP 5 : Wait Two Minutes

Wait two minutes for the results. If the colour change extends beyond the red line on the tube, you're over the limit.



COMPANY PROFILE



Bellerophon Biotec Pvt. Ltd. was founded in 2012.

Company Philosophy: To provide high quality diagnostic reagents, instruments and systems with technical support to the Health Services of India.

Our Commitment: To provide the finest quality products, matched with friendly and efficient after-sales support.

Technical Assistance: Experienced sales personnel are always available for direction and guidance. In-house technical support, both before and after the sale is just a phone call away.

We strive to make the Company customer-orientated by focusing on our customers' needs, ensuring that the products we distribute and services we provide are of the highest quality.

SensiDrive-

A Disposable Breath Analyzer— is a unique product launched by us, which is poised to drastically reduce road accidents in India, thus saving scores of precious lives.

Regd. Office : 57-A (Basement), Sant Nagar, East of Kailash, New Delhi-110065 (India)

Address for Communication: 34, Knowledge Park-1, Greater Noida (UP) 201310

Tel.: +91-120-2326042, +91-11-40701328 Fax: 0120-2326078

E-mail : info@bellerophonbiotec.com

Website : www.bellerophonbiotec.com (A Chimera group of Company)